



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Dog Serum NMB-NDS

Article Name	Dog Serum
Biozol Catalog Number	NMB-NDS
Supplier Catalog Number	NDS
Alternative Catalog Number	NMB-NDS
Manufacturer	NordicMubio
Category	Zellen/Zellkultur
Species Reactivity	Canine
Product Description	As a qualitative serum with an optimized protein content to study the normal protein composition of dog serum. To identify or quantitate a normal serum protein component using a variety of immunodiffusion techniques including immunoelectrophoresis, s
Buffer	Delipidated, lyophilized, stable reference serum. No preservative added.
Source	Pooled normal Dog Serum. The NDS dog serum was obtained from normal, healthy dogs of the Bouvier breed. Regulations are in place in the Netherlands for these dogs to have had some form of vaccination. The Normal Dog Serum (NDS) was derived from dogs that
Application Notes	Normal serum protein reference sera are shipped lyophilized at ambient temperature. After arrival store at +4C. Prolonged storage may be at -20C. After opening the vial, the lyophilized content is reconstituted by adding 1 ml of sterile distilled water