

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

FIX&PERM Kit 1000 Cell Fixation and Permeabilization Kit (CE/IVD) NMB-GAS-002-1

Article Name	FIX&PERM Kit 1000 Cell Fixation and Permeabilization Kit (CE/IVD)
Biozol Catalog Number	NMB-GAS-002-1
Supplier Catalog Number	GAS-002-1
Alternative Catalog Number	NMB-GAS-002-1
Manufacturer	NordicMubio
Category	Sonstiges
Application	FC
Product Description	Flow cytometric analyses with monoclonal antibodies were so far mainly restricted to cell surface molecules. Intracellular structures such as cytoplasmic or nuclear enzymes, oncoproteins, cytokines, immunoglobulins etc. were largely excluded from suc...
Buffer	This FIX&PERM Cell Fixation and Permeabilization Kit contains 2 reagents: Fixation Medium (Reagent A) and Permeabilization Medium (Reagent B). It is intended for first fixing cells in suspension with Reagent A and then permeabilizing the cell membrane
Application Notes	Biological fluids (blood, bone marrow, and others) must be collected under sterile conditions. Anticoagulation with EDTA or heparin is recommended. The samples should be stored at room temperature until used. For optimal results, samples should be proces