

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### FIX&PERM Sample Kit (CE/IVD) NMB-GAS-002M

Article Name	FIX&PERM Sample Kit (CE/IVD)
Biozol Catalog Number	NMB-GAS-002M
Supplier Catalog Number	GAS-002M
Alternative Catalog Number	NMB-GAS-002M
Manufacturer	NordicMubio
Category	Sonstiges
Application	FC
Product Description	Flow cytometric analyses with monoclonal antibodies were so far mainly restricted to cell surface molecules. Intracellular structures such as cytoplasmic or nuclear enzymes, oncoproteins, cytokines, immunoglobulins etc. were largely excluded from suc...
Buffer	This FIX&PERM Cell Fixation and Permeabilization Kit contains 2 reagents: Fixation Medium (Reagent A) and Permeabilization Medium (Reagent B). It is intended for first fixing cells in suspension with Reagent A and then permeabilizing the cell membrane
Application Notes	Biological fluids (blood, bone marrow, and others) must be collected under sterile conditions. Anticoagulation with EDTA or heparin is recommended. The samples should be stored at room temperature until used. For optimal results, samples should be proces