

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Goat anti Dog Serum proteins, Clone: [Polyclonal], Monoclonal NMB-GAD/TSP

Article Name	Goat anti Dog Serum proteins, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-GAD/TSP
Supplier Catalog Number	GAD/TSP
Alternative Catalog Number	NMB-GAD/TSP
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Product Description	In immunoelectrophoresis against pooled serum and enriched serum proteins fractions precipitation can be observed of not less than 20 individual proteins components. In precipitating techniques as immunoelectrophoresis and radial immunodiffusion (Ouc...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added. Reconstitute the lyophilized
Source	Pooled whole dog serum and partly purified serum fractions. Freund's complete adjuvant is used in the first step of the immunization procedure.

Formula	Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added.
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 $\mu$ l serum or equivalent against 120 $\mu$ l antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 $\mu$ l antiserum in 3 mm diameter center well and 2 $\mu$ l serum samples (neat