

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Human IgG IgA IgM IgD (heavy and light chains), Clone: [Polyclonal], Monoclonal NMB-GAHU/IG(GAMD)/7S

| Article Name | Goat anti Human IgG IgA IgM IgD (heavy and light chains), Clone: [Polyclonal], Monoclonal |
|----------------------------|---|
| Biozol Catalog Number | NMB-GAHU/IG(GAMD)/7S |
| Supplier Catalog Number | GAHu/Ig(GAMD)/7S |
| Alternative Catalog Number | NMB-GAHU/IG(GAMD)/7S |
| Manufacturer | NordicMubio |
| Host | Goat |
| Category | Antikörper |
| Product Description | The defined antibody specificity is directed to the major isotypes of the human immunoglobulin system (classes and both light chain types) including antibody to common determinants and to class specific determinants and to the surface determinants of |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove i |
| Source | Highly purified IgG, and pools of homogenous IgA, IgM and IgD isolated from human serum. Freunds complete adjuvant is used in the first step of the immunization procedure. |

| | Formula | Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. | |
|--|-------------------|--|--|
| | Application Notes | Indirect immunofluorescence, ELISA, Dot blot, Immunoblotting. | |