

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Monkey IgM (Fc specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal NMB-GAMON/IGM(FC)/TRITC

| Article Name | Goat anti Monkey IgM (Fc specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-GAMON/IGM(FC)/TRITC |
| Supplier Catalog Number | GAMon/IgM(Fc)/TRITC |
| Alternative Catalog Number | NMB-GAMON/IGM(FC)/TRITC |
| Manufacturer | NordicMubio |
| Host | Goat |
| Category | Antikörper |
| Conjugation | TRITC |
| Product Description | The reactivity of the antiserum is directed to the Fc subunit of the IgM molecule which expresses strict isotypic (class) specificity. It does not react with any non-Ig protein in monkey serum, as tested by immunoelectrophoresis and double radial imm |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun do |

| Source | Purified polyclonal IgM and pools of homogenous IgM isolated from pooled monkey serum. Freunds complete adjuvant is used in the first step of the immunization procedure. |
|-------------------|--|
| Formula | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Application Notes | ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence. |