

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Goat anti Mouse IgG2 (subclass specific), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-GAM/IGG2/PO

Article Name	Goat anti Mouse IgG2 (subclass specific), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal
Biozol Catalog Number	NMB-GAM/IGG2/PO
Supplier Catalog Number	GAM/IgG2/PO
Alternative Catalog Number	NMB-GAM/IGG2/PO
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Conjugation	HRP
Product Description	The reactivity of the antiserum is directed to the subclass IgG2a and IgG2b. It reacts with determinants specific for IgG2a and IgG2b and iso and allospecific determinants. It does not react with other subclasses of IgG, IgG/Fab fragments, IgM and I
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, sp

Source	Pools of purified homogenous IgG2a and IgG2b isolated from pooled mouse sera belonging to the allelic types Igh-1a and Igh-1b. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.