

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Mouse IgG2 (subclass specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal NMB-GAM/IGG2/TRITC

| Article Name | Goat anti Mouse IgG2 (subclass specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-GAM/IGG2/TRITC |
| Supplier Catalog Number | GAM/IgG2/TRITC |
| Alternative Catalog Number | NMB-GAM/IGG2/TRITC |
| Manufacturer | NordicMubio |
| Host | Goat |
| Category | Antikörper |
| Conjugation | TRITC |
| Product Description | The reactivity of the antiserum is directed to both subclasses IgG2a and IgG2b. It reacts with determinants specific for IgG2a and IgG2b and iso and allospecific determinants. It does not react with other subclasses of IgG, IgG/Fab fragments, IgM and |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun do |

| Source | Pools of purified homogenous IgG2a and IgG2b isolated from mouse serum of strains belonging to the allelic types Igh-1a and Igh-1b. Freunds complete adjuvant is used in the first step of the immunization procedure. |
|-------------------|--|
| Formula | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Application Notes | ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence. |