

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Rabbit Serum proteins , conjugated with FITC, Clone: [Polyclonal], Monoclonal NMB-GAR/TSP/FITC

Article Name	Goat anti Rabbit Serum proteins , conjugated with FITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-GAR/TSP/FITC
Supplier Catalog Number	GAR/TSP/FITC
Alternative Catalog Number	NMB-GAR/TSP/FITC
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Conjugation	FITC
Product Description	Tested in immunoelectrophoresis and double radial immunodiffusion against pooled serum precipitation of at least 14 individual proteins have been observed. Direct fluorescent staining of rabbit serum proteins at the cellular and subcellular level of
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	FITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into sm

Source	Pooled whole serum and selected protein fractions isolated from pooled rabbit serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	FITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.