

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Rat IgE (Fc specific), Clone: [Polyclonal], Monoclonal NMB-GARA/IGE(FC)/7S

Goat anti Rat IgE (Fc specific), Clone: [Polyclonal], Monoclonal
NMB-GARA/IGE(FC)/7S
GARa/IgE(Fc)/7S
NMB-GARA/IGE(FC)/7S
NordicMubio
Goat
Antikörper
The reactivity of the antiserum is directed to the Fc subunits of the IgE molecule which represents strict isotype (class) specificity. In immunoelectrophoresis and double radial immunodiffusion, using various antiserum concentrations against appropr
Monoclonal
[Polyclonal]
Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. Reconstitute the lyophilized antiserum by adding 1 ml sterile distilled water.
Purified homogenous IgE isolated from rat serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.

Indirect immunofluorescence, ELISA, Dot blot, Immunoblotting.