

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Rat IgG2a (subclass specific), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-GARA/IGG2A/PO

| | |
|----------------------------|---|
| Article Name | Goat anti Rat IgG2a (subclass specific), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal |
| Biozol Catalog Number | NMB-GARA/IGG2A/PO |
| Supplier Catalog Number | GARa/IgG2a/PO |
| Alternative Catalog Number | NMB-GARA/IGG2A/PO |
| Manufacturer | NordicMubio |
| Host | Goat |
| Category | Antikörper |
| Species Reactivity | Rat |
| Conjugation | HRP |
| Product Description | The reactivity of the antiserum is directed to the subclass IgG2a. It does not react with other subclasses of IgG, IgG/Fab fragments, IgM and IgA or any non-Ig protein in rat serum, as tested by immunoelectrophoresis and double radial immunodiffusion... |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. Reconstitute the lyophilized antiserum by adding 1 ml sterile |

| | |
|-------------------|---|
| Source | Pools of purified homogenous IgG2a isolated from pooled rat serum. Freund's complete adjuvant is used in the first step of the immunization procedure. |
| Formula | Peroxidase-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Application Notes | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting. |