

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Aspergillus niger Amyloglucosidase, Clone: [Polyclonal], Monoclonal NMB-NE018/PAB

| Article Name               | Rabbit anti Aspergillus niger Amyloglucosidase, Clone: [Polyclonal],<br>Monoclonal   |
|----------------------------|--|
| Biozol Catalog Number      | NMB-NE018/PAB  |
| Supplier Catalog Number    | NE018/PAb  |
| Alternative Catalog Number | NMB-NE018/PAB  |
| Manufacturer               | NordicMubio  |
| Host                       | Rabbit   |
| Category                   | Antikörper   |
| Product Description        | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and     |
| Clonality                  | Monoclonal   |
| Clone Designation          | [Polyclonal]   |
| Buffer                     | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut |
| Source                     | Amyloglucosidase isolated and purified from Aspergillus niger.<br>Freunds complete adjuvant is used in the first step of the<br>immunization procedure.  |

| Formula           | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
|-------------------|--|
| Application Notes | Immunoprecipitation.   |