

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Bacillus polymyxa Protease neutral, Clone: [Polyclonal], Monoclonal NMB-NE143/PAB

| Article Name | Rabbit anti Bacillus polymyxa Protease neutral, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-NE143/PAB |
| Supplier Catalog Number | NE143/PAb |
| Alternative Catalog Number | NMB-NE143/PAB |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Product Description | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| UniProt | P29148 |
| Buffer | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut |

| Source | Neutral protease isolated and purified from Bacillus polymyxa. Freunds complete adjuvant is used in the first step of the immunization procedure. |
|-------------------|--|
| Formula | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Immunoprecipitation. |