

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Bacillus subtilis a-Amylase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE015/BIO

Article Name	Rabbit anti Bacillus subtilis a-Amylase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE015/BIO
Supplier Catalog Number	NE015/Bio
Alternative Catalog Number	NMB-NE015/BIO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Application	DOT, ELISA, ICC, WB
Species Reactivity	Bacteria
Conjugation	Biotin
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	P00691

Buffer	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile di
Source	alpha-Amylase isolated and purified from Bacillus subtilis. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.