

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti Bovine albumin, conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-RAB/ALB/PO**

|                            |   |
|----------------------------|---|
| Article Name               | Rabbit anti Bovine albumin, conjugated with Horseradish peroxidase,<br>Clone: [Polyclonal], HRP, Monoclonal   |
| Biozol Catalog Number      | NMB-RAB/ALB/PO  |
| Supplier Catalog Number    | RAB/Alb/PO  |
| Alternative Catalog Number | NMB-RAB/ALB/PO  |
| Manufacturer               | NordicMubio   |
| Host                       | Rabbit  |
| Category                   | Antikörper  |
| Conjugation                | HRP   |
| Product Description        | Tested in immunoelectrophoresis and double radial immunodiffusion against pooled normal bovine serum and purified bovine albumin. One characteristic precipitin line is obtained against pooled normal bovine serum using different antigen/antibody conce... |
| Clonality                  | Monoclonal  |
| Clone Designation          | [Polyclonal]  |
| Buffer                     | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, |

|                   |   |
|-------------------|---|
| Source            | Purified albumin isolated from pooled bovine serum. Freund's complete adjuvant is used in the first step of the immunization procedure.   |
| Formula           | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Application Notes | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.  |