

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Bovine IgM (Fc specific), conjugated with FITC, Clone: [Polyclonal], Monoclonal NMB-RAB/IGM(FC)/FITC

Article Name	Rabbit anti Bovine IgM (Fc specific), conjugated with FITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RAB/IGM(FC)/FITC
Supplier Catalog Number	RAB/IgM(Fc)/FITC
Alternative Catalog Number	NMB-RAB/IGM(FC)/FITC
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	FITC
Product Description	The reactivity of the antiserum is directed to the Fc subunit of the IgM molecule which expresses strict isotypic (class) specificity. It does not react with any non-Ig protein in bovine serum, as tested by immunoelectrophoresis and double radial imm
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Fluorescein-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water,

Source	Purified normal IgM isolated from pooled bovine serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Fluorescein-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.