

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Bovine Uridine-5-diphosphoglucose dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE189/BIO

| Article Name | Rabbit anti Bovine Uridine-5-diphosphoglucose dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-NE189/BIO |
| Supplier Catalog Number | NE189/Bio |
| Alternative Catalog Number | NMB-NE189/BIO |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Conjugation | Biotin |
| Product Description | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| UniProt | P12379 |
| Buffer | Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile di |

| Source | Uridine-5-diphosphoglucosedehydrogenase isolated and purified from bovine liver. Freunds complete adjuvant is used in the first step of the immunization procedure. |
|-------------------|--|
| Formula | Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting. |