



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Bovine Uridine-5-diphosphoglucose pyrophosphorylase, Clone: [Polyclonal], Monoclonal NMB-NE190/7S

| Article Name | Rabbit anti Bovine Uridine-5-diphosphoglucose pyrophosphorylase, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-NE190/7S |
| Supplier Catalog Number | NE190/7S |
| Alternative Catalog Number | NMB-NE190/7S |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Product Description | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| UniProt | Q07130 |
| Buffer | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distille |

| Source | Uridine-5-diphosphoglucose pyrophosphorylase isolated and purified from bovine liver. Freunds complete adjuvant is used in the first step of the immunization procedure. |
|-------------------|--|
| Formula | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting. |