

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Crotalus durissus terrificus venom Phosphodiesterase I, Clone: [Polyclonal], Monoclonal NMB-NE126/PAB

Article Name	Rabbit anti Crotalus durissus terrificus venom Phosphodiesterase I, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE126/PAB
Supplier Catalog Number	NE126/PAb
Alternative Catalog Number	NMB-NE126/PAB
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut
Source	Phosphodiesterase isolated and purified from Crotalus durissus terrificus venom. Freunds complete adjuvant is used in the first step of the immunization procedure.

Formula	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	Immunoprecipitation.