

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Dog Serum proteins, conjugated with FITC, Clone: [Polyclonal], Monoclonal NMB-RAD/TSP/FITC

Article Name	Rabbit anti Dog Serum proteins, conjugated with FITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RAD/TSP/FITC
Supplier Catalog Number	RAD/TSP/FITC
Alternative Catalog Number	NMB-RAD/TSP/FITC
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	FITC
Product Description	Tested in immunoelectrophoresis and double radial immunodiffusion against pooled serum precipitation of at least 12 individual proteins have been observed. Direct fluorescent staining of dog serum proteins at the cellular and subcellular level of app
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	FITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into

Source	Pooled whole serum and selected protein fractions isolated from pooled dog serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	FITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.