

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Enterobacter aerogenes Glycerol dehydrogenase, Clone: [Polyclonal], Monoclonal NMB-NE068/PAB

Article Name	Rabbit anti Enterobacter aerogenes Glycerol dehydrogenase, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE068/PAB
Supplier Catalog Number	NE068/PAb
Alternative Catalog Number	NMB-NE068/PAB
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	G0EAP7
Buffer	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut

Source	Glycerol dehydrogenase isolated and purified from Enterobacter aerogenes. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	Immunoprecipitation.