

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Human fibonectin, Clone: [Polyclonal], Monoclonal NMB-RAHU/FBN

Article Name	Rabbit anti Human fibonectin, Clone: [Polyclonal], Monoclonal
Audie Haille	Nappic and Human inconectin, Cione. [Folyclonal], Monoclonal
Biozol Catalog Number	NMB-RAHU/FBN
Supplier Catalog Number	RAHu/Fbn
Alternative Catalog Number	NMB-RAHU/FBN
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	The reactivity of the antiserum is restricted to fibronectin. In immunoelectrophoresis and radial immunodiffusion, using various antiserum concentrations against human plasma, sera and purified fibronectin a single precipitin line is obtained. Testin
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	P02751
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum. Sodium azide 0.7 mg/ml Total protein and IgG concentrations in the antiserum are comparable to those of normal pooled rabbit serum. No foreign proteins added. It is reconstituted by addi
Source	Fibrinogen (clotting factor I) is a heat labile beta glycoprotein (molecular weight 340,000) and consists of three pairs of chains bound by disulphide bonds. It is synthesized in hepatocytes under genetic control. It is the precursor of fibrin, which is

Formula	Delipidated, heat inactivated, lyophilized, stable whole antiserum. Sodium azide 0.7 mg/ml Total protein and lgG concentrations in the antiserum are comparable to those of normal pooled rabbit serum. No foreign proteins added.
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 μ l amniotic fluid or equivalent against 120 μ l antiserum. In double radial immunodiffusion use a rosette arrangement with 10 μ l antiserum in 3 mm diameter center well and 2 μ l amniotic fluid or hepatom