

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Mouse IgG Fab, conjugated with TRITC, Clone: [Polyclonal], Monoclonal NMB-RAM/FAB/TRITC

Article Name	Rabbit anti Mouse IgG Fab, conjugated with TRITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RAM/FAB/TRITC
Supplier Catalog Number	RAM/Fab/TRITC
Alternative Catalog Number	NMB-RAM/FAB/TRITC
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	TRITC
Product Description	The reactivity of the antiserum is directed to the Fab subunits of intact IgG, IgA IgM and other Ig classes of both light chain types, with their Fab or F(ab)2 subunits, and with free light chains of kappa and lambda type. In immunoelectrophoresis an
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	TRITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun
Source	Purified Fab of normal mouse IgG. Freunds complete adjuvant is used in the first step of the immunization procedure.

Formula	TRITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.