

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti Mouse Serum proteins (non-Ig), Clone: [Polyclonal], Monoclonal NMB-RAM/AGGS**

Article Name	Rabbit anti Mouse Serum proteins (non-Ig), Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RAM/AGGS
Supplier Catalog Number	RAM/AGGS
Alternative Catalog Number	NMB-RAM/AGGS
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	In immunoelectrophoresis against pooled serum and enriched serum proteins fractions precipitation can be observed of not less than 15 individual proteins components. In precipitating techniques as immunoelectrophoresis and radial immunodiffusion (Ouc...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added. Reconstitute the lyophilized
Source	Pooled immunoglobulin-depleted mouse serum and partly purified serum fractions. Freund's complete adjuvant is used in the first step of the immunization procedure.

Formula	Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added.
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat