

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti pigeon IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Monoclonal NMB-RAP/IGG(H+L)/FITC**

Article Name	Rabbit anti pigeon IgG (heavy and light chains), conjugated with FITC, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RAP/IGG(H+L)/FITC
Supplier Catalog Number	RAP/IgG(H+L)/FITC
Alternative Catalog Number	NMB-RAP/IGG(H+L)/FITC
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	FITC
Product Description	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in pigeon serum, a...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	FITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added It is reconstituted by adding 2 ml sterile distilled water, spun down to remove insoluble particles, divided into s

Source	Purified normal IgG isolated from pooled pigeon serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formula	FITC-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added
Application Notes	ELISA, Immunocytochemistry, Immunohistochemistry (frozen), (In)direct immunofluorescence.