

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Porcine Amino acid arylamidase, Clone: [Polyclonal], Monoclonal NMB-NE011/PAB

| Article Name | Rabbit anti Porcine Amino acid arylamidase, Clone: [Polyclonal], Monoclonal |
|----------------------------|---|
| Biozol Catalog Number | NMB-NE011/PAB |
| Supplier Catalog Number | NE011/PAb |
| Alternative Catalog Number | NMB-NE011/PAB |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Product Description | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut |
| Source | Amino acid arylamidase isolated and purified from porcine kidney. Freunds complete adjuvant is used in the first step of the immunization procedure. |

| Formula | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
|-------------------|--|
| Application Notes | Immunoprecipitation. |