

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Pseudomonas fluorescens b-Galactosidase dehydrogenase, Clone: [Polyclonal], Monoclonal NMB-NE170/7S

Article Name	Rabbit anti Pseudomonas fluorescens b-Galactosidase dehydrogenase, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE170/7S
Supplier Catalog Number	NE170/7S
Alternative Catalog Number	NMB-NE170/7S
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	P11886
Buffer	Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distille

Source	beta-Galactose dehydrogenase is isolated and purified from Pseudomonas fluorescens. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting.