

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Rat transferrin, Clone: [Polyclonal], Monoclonal NMB-RARA/TRF

Article Name	Rabbit anti Rat transferrin, Clone: [Polyclonal], Monoclonal	
Biozol Catalog Number	NMB-RARA/TRF	
Supplier Catalog Number	RARa/Trf	
Alternative Catalog Number	NMB-RARA/TRF	
Manufacturer	NordicMubio	
Host	Rabbit	
Category	Antikörper	
Product Description	The reactivity of the antiserum is restricted to transferrin as tested in immunoelectrophoresis and radial immunodiffusion. A characteristic single precipitin line is obtained with normal serum which shows a reaction of identity with the precipitin l	
Clonality	Monoclonal	
Clone Designation	[Polyclonal]	
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added. Reconstitute the lyophilize	
Source	Purified transferrin isolated from pooled rat serum. Freunds complete adjuvant is used in the first step of the immunization procedure.	
Formula	Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added.	

Annl	ication	Notes
ADDI	icalion	MOLES

Precipitation assays. In immunoelectrophoresis use 2 μl serum or equivalent against 120 μl antiserum. In double radial immunodiffusion use a rosette arrangement with 10 μl antiserum in 3 mm diameter centre well and 2 μl serum samples (neat and serially d