

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Rhodopseudomonas spheroides b-Hydroxybutyrate dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE078/BIO

Article Name	Rabbit anti Rhodopseudomonas spheroides b-Hydroxybutyrate dehydrogenase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE078/BIO
Supplier Catalog Number	NE078/Bio
Alternative Catalog Number	NMB-NE078/BIO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	Biotin
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile di

Source	beta-Hydroxybutyrate dehydrogenase isolated and purified from Rhodopseudomonas spheroides. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.