

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Saccharomyces cerevisiae 3-Phosphoglyceric phosphokinase, Clone: [Polyclonal], Monoclonal NMB-NE130/PAB

Article Name	Rabbit anti Saccharomyces cerevisiae 3-Phosphoglyceric phosphokinase, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-NE130/PAB
Supplier Catalog Number	NE130/PAb
Alternative Catalog Number	NMB-NE130/PAB
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Product Description	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut
Source	3-Phosphoglyceric phosphokinase isolated and purified from bakers yeast. Freunds complete adjuvant is used in the first step of the immunization procedure

Formula	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Application Notes	Immunoprecipitation.