

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Sheep IgM (heavy and light chains), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-RASH/IGM(H+L)/PO

Article Name	Rabbit anti Sheep IgM (heavy and light chains), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal
Biozol Catalog Number	NMB-RASH/IGM(H+L)/PO
Supplier Catalog Number	RASh/IgM(H+L)/PO
Alternative Catalog Number	NMB-RASH/IGM(H+L)/PO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	HRP
Product Description	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgM molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in sheep serum, as
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Horseradish peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particl

Source	Purified normal IgM isolated from pooled sheep serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Horseradish peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.