

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Sheep Sorbitol dehydrogenase, Clone: [Polyclonal], Monoclonal NMB-NE149/PAB

| | |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name | Rabbit anti Sheep Sorbitol dehydrogenase, Clone: [Polyclonal], Monoclonal |
| Biozol Catalog Number | NMB-NE149/PAB |
| Supplier Catalog Number | NE149/PAb |
| Alternative Catalog Number | NMB-NE149/PAB |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Product Description | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and... |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| UniProt | P07846 |
| Buffer | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted |

| | |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Source | Sorbitol dehydrogenase isolated and purified from sheep liver. Freund's complete adjuvant is used in the first step of the immunization procedure. |
| Formula | Purified hyperimmune rabbit IgG antibodies lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Immunoprecipitation. |