

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit anti Swine transferrin, Clone: [Polyclonal], Monoclonal NMB-RASW/TRF/7S

| Article Name | Rabbit anti Swine transferrin, Clone: [Polyclonal], Monoclonal |
|----------------------------|--|
| Biozol Catalog Number | NMB-RASW/TRF/7S |
| Supplier Catalog Number | RASw/Trf/7S |
| Alternative Catalog Number | NMB-RASW/TRF/7S |
| Manufacturer | NordicMubio |
| Host | Rabbit |
| Category | Antikörper |
| Product Description | The reactivity of the antiserum is restricted to transferrin as tested in immunoelectrophoresis and radial immunodiffusion. A characteristic single precipitin line is obtained with normal serum which shows a reaction of identity with the precipitin l |
| Clonality | Monoclonal |
| Clone Designation | [Polyclonal] |
| Buffer | Purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. Reconstitute the lyophilized antiserum by adding 1 ml sterile distilled water. |
| Source | Purified transferrin isolated from pooled swine serum. Freunds complete adjuvant is used in the first step of the immunization procedure. |
| Formula | Purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. |

Indirect immunofluorescence, ELISA, Dot blot, Immunoblotting.