

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rabbit anti Swine transferrin, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-RASW/TRF/BIO

Article Name	Rabbit anti Swine transferrin, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-RASW/TRF/BIO
Supplier Catalog Number	RASw/Trf/Bio
Alternative Catalog Number	NMB-RASW/TRF/BIO
Manufacturer	NordicMubio
Host	Rabbit
Category	Antikörper
Conjugation	Biotin
Product Description	The reactivity of the antiserum is restricted to transferrin as tested in immunoelectrophoresis and radial immunodiffusion. A characteristic single precipitin line is obtained with normal serum which shows a reaction of identity with the precipitin l
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Biotin-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added Reconstitute the lyophilized antiserum by adding 1 ml sterile distilled water.

Source	Purified transferrin isolated from pooled swine serum. Freunds complete adjuvant is used in the first step of the immunization procedure.
Formula	Biotin-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting.