

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

HPV-16 E6 + HPV-18 E6 Antibody, Clone: [HPV16/1295 + HPV18/1297], Mouse, Monoclonal NSJ-V8334SAF-100UG

| | |
|----------------------------|---|
| Article Name | HPV-16 E6 + HPV-18 E6 Antibody, Clone: [HPV16/1295 + HPV18/1297], Mouse, Monoclonal |
| Biozol Catalog Number | NSJ-V8334SAF-100UG |
| Supplier Catalog Number | V8334SAF-100UG |
| Alternative Catalog Number | NSJ-V8334SAF-100UG |
| Manufacturer | NSJ Bioreagents |
| Host | Mouse |
| Category | Antikörper |
| Application | ELISA |
| Species Reactivity | Virus |
| Immunogen | HPV18 E6-i-galactosidase fusion protein was used as the immunogen for this HPV-16 E6 + HPV-18 E6 antibody cocktail. |
| Product Description | Human papilloma viruses (HPVs) can be classified as either high risk or low risk according to their association with cancer. HPV16 and HPV18 are the most common of the high risk group while HPV6 and HPV11 are among the low risk types. Approximately 9... |
| Clonality | Monoclonal |
| Concentration | 1 mg/ml in 1X PBS, BSA free, sodium azide free |
| Clone Designation | [HPV16/1295 + HPV18/1297] |
| Isotype | Mouse IgG1, kappa |

| | |
|--------------------|--|
| Buffer | 1 mg/ml in 1X PBS, BSA free, sodium azide free |
| Purity | Protein G affinity chromatography |
| Form | 1 mg/ml in 1X PBS, BSA free, sodium azide free |
| Target | HPV-16 E6 + HPV-18 E6 |
| Antibody Type | Primary Antibody |
| Application Dilute | ELISA (order BSA-free format for coating) |
| Application Notes | Optimal dilution of the antibody should be determined by the researcher. |