



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Ab SpinTrap columns CRM-AB-001P-50

Article Name	Ab SpinTrap columns
Biozol Catalog Number	CRM-AB-001P-50
Supplier Catalog Number	Ab-001P
Alternative Catalog Number	CRM-AB-001P-50
Manufacturer	Creative BioMart
Category	Sonstiges
Product Description	Columns can be used with a standard microcentrifuge and one purification run takes less than 20 min Designed for efficient, small-scale purification of monoclonal and polyclonal antibodies from serum and cell culture supernatants
Application Notes	Columns can be used with a standard microcentrifuge and one purification run takes less than 20 min Designed for efficient, small-scale purification of monoclonal and polyclonal antibodies from serum and cell culture supernatants.