

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Streptavidin Sepharose High Performance affinity resin CRM-ST-383C-5

Article Name	Streptavidin Sepharose High Performance affinity resin
Biozol Catalog Number	CRM-ST-383C-5
Supplier Catalog Number	St-383C
Alternative Catalog Number	CRM-ST-383C-5
Manufacturer	Creative BioMart
Category	Sonstiges
Product Description	The interaction between streptavidin and biotin is very strong and requires denaturing conditions for elution. This strong interaction can be utilized in the purification of antigens. Biotinylated antibody-antigen complexes bind to the HiTrap Strepta...
Application Notes	The interaction between streptavidin and biotin is very strong and requires denaturing conditions for elution. This strong interaction can be utilized in the purification of antigens. Biotinylated antibody-antigen complexes bind to the HiTrap Streptavidi