

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

DFF40/CAD polyclonal antibody ENZ-ADI-AAP-350-E

| DFF40/CAD polyclonal antibody |
|---|
| ENZ-ADI-AAP-350-E |
| ADI-AAP-350-E |
| ENZ-ADI-AAP-350-E-100 |
| Enzo Life Sciences |
| Rabbit |
| Antikörper |
| IHC, WB |
| Mouse, Rat |
| Synthetic peptide corresponding to the sequence near the middle of mouse DFF40/CAD. |
| DFF40/CAD polyclonal antibody |
| Polyclonal |
| O54788 |
| Peptide affinity purified. |
| Liquid. In PBS containing 0.02% sodium azide. |
| DFF |
| Detects a band of ~43kDa by Western blot. |
| |