

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Ubiquitin, (agarose immobilized) ENZ-BML-UW8630-0500

Article Name	Ubiquitin, (agarose immobilized)
Biozol Catalog Number	ENZ-BML-UW8630-0500
Supplier Catalog Number	BML-UW8630-0500
Alternative Catalog Number	ENZ-BML-UW8630-0500-0.5
Manufacturer	Enzo Life Sciences
Category	Zellen/Zellkultur
Product Description	Ubiquitin, (agarose immobilized)...
UniProt	P0CG47
Source	Ubiquitin is coupled to agarose, via the epsilon amino groups of its internal lysine residues, leaving the C-terminus free. The loading is approximately 10mg protein/ml of settled resin.
Form	Liquid. Suspension.
Target	Ubiquitin
Application Notes	Ubiquitin-agarose can be used to affinity purify ubiquitin-binding proteins.