

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Cytokeratin 18 (KRT18) Antibody, Clone: [KRT18/7056R], Rabbit, Monoclonal NBT-3875-RBM18-P0

| Article Name               | Recombinant Cytokeratin 18 (KRT18) Antibody, Clone: [KRT18/7056R], Rabbit, Monoclonal  |
|----------------------------|--|
| Biozol Catalog Number      | NBT-3875-RBM18-P0  |
| Supplier Catalog Number    | 3875-RBM18-P0  |
| Alternative Catalog Number | NBT-3875-RBM18-P0-20,NBT-3875-RBM18-P0-100   |
| Manufacturer               | NeoBiotechnologies   |
| Host                       | Rabbit   |
| Category                   | Antikörper   |
| Application                | IHC  |
| Species Reactivity         | Human  |
| Immunogen                  | Recombinant human full-length KRT18 protein  |
| Product Description        | This MAb reacts with a wide variety of simple epithelia. It does not react with stratified squamous epithelia. It reacts with epithelial tumors of the gastrointestinal tract, lung, breast, pancreas, ovary, and thyroid. Cytokeratin 18, which belongs t |
| Clonality                  | Monoclonal   |
| Clone Designation          | [KRT18/7056R]  |
| Molecular Weight           | 45kDa  |
| NCBI                       | 3875   |

| UniProt           | P05783   |
|-------------------|--|
| Form              | 200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are a side. Also available WITHOUT BSA & Double at 1.0mg/ml.  |
| Antibody Type     | Recombinant Monoclonal Antibody  |
| Application Notes | Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95&degC followed by cooling at RT for 20 minutes),Optimal dilutio |