

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) Antibody, IgG1, Clone: [SPM168], Mouse, Monoclonal NBT-800-MSM1X-P0

Article Name	Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) Antibody, IgG1, Clone: [SPM168], Mouse, Monoclonal
Biozol Catalog Number	NBT-800-MSM1X-P0
Supplier Catalog Number	800-MSM1X-P0
Alternative Catalog Number	NBT-800-MSM1X-P0-20,NBT-800-MSM1X-P0-100
Manufacturer	NeoBiotechnologies
Host	Mouse
Category	Antikörper
Application	FC, IF, IHC, WB
Species Reactivity	Human
Immunogen	Crude human uterus extract
Product Description	Recognizes a protein of 150kDa, which is identified as the high molecular weight variant of Caldesmon. Two closely related variants of human caldesmon have been identified which are different in their electrophoretic mobility and cellular distributio
Clonality	Monoclonal
Clone Designation	[SPM168]
Molecular Weight	150kDa
Isotype	lgG1

NCBI	800
UniProt	Q05682
Form	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are also available WITHOUT BSA & Double at 1.0mg/ml.
Antibody Type	Monoclonal Antibody
Application Notes	Western Blot (2-4ug/ml), Flow Cytometry (1-2ug/million cells), Immunofluorescence (1-2ug/ml), Immunohistochemistry (Formalinfixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 1mM EDTA, pH 7.5-8.