

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Penicillin (Group of beta-Lactam Antibiotics) Antibody, IgG1, Clone: [Pen-9], Mouse, Monoclonal NBT-MSM1-231-P1ABX

Article Name	Penicillin (Group of beta-Lactam Antibiotics) Antibody, IgG1, Clone: [Pen-9], Mouse, Monoclonal
Biozol Catalog Number	NBT-MSM1-231-P1ABX
Supplier Catalog Number	MSM1-231-P1ABX
Alternative Catalog Number	NBT-MSM1-231-P1ABX-100
Manufacturer	NeoBiotechnologies
Host	Mouse
Category	Antikörper
Application	ELISA
Immunogen	Penicilloyl-transferrin conjugate
Product Description	This antibody reacts mainly with the hiazolidine ring of penicillin, but not with the lactam ring. May react with the fused beta lactam/thiazolidine ring of the penicilloyl group and would appear sensitive to modification of the shared N atom of the
Clonality	Monoclonal
Clone Designation	[Pen-9]
Molecular Weight	356 Daltons
Isotype	lgG1
NCBI	Not Applicable

Form	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & Double are applied to the substitution of Ab purified from Bioreactor Concentrate by Protein A/G. Protei
Antibody Type	Monoclonal Antibody
Application Notes	ELISA (For coating, use Ab at 1-5ug/ml, order Ab without BSA)Optimal dilution for a specific application should be determined.