



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Triple FLAG Peptide, Unconjugated ROC-000-000-H47

Article Name	Triple FLAG Peptide, Unconjugated
Biozol Catalog Number	ROC-000-000-H47
Supplier Catalog Number	000-000-Н47
Alternative Catalog Number	ROC-000-000-H47
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE
Conjugation	Unconjugated
Concentration	1.0 mg/ml by dry weight
Buffer	None
Form	Lyophilized
Application Dilute	IP: YES - see PROTOCOL, WB: YES as blocking peptide in PCA
Application Notes	The Triple FLAG Peptide has been tested by SDS-PAGE. A working concentration of 500 μ g/mL for triple FLAG peptide is typically used for competitive elution of FLAG containing proteins from anti-DYKDDDDK affinity gels (p/n 200-350-383). For use as a co