

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Surface Lipoprotein p27 Control Protein, Unconjugated ROC-000-001-C30

Article Name	Surface Lipoprotein p27 Control Protein, Unconjugated
Biozol Catalog Number	ROC-000-001-C30
Supplier Catalog Number	000-001-C30
Alternative Catalog Number	ROC-000-001-C30
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE, WB
Conjugation	Unconjugated
Concentration	1.0 mg/mL by modified Lowry assay
NCBI	1194336
UniProt	O50951
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Source	Borrelia burgdorferi
Form	Liquid (sterile filtered)
Application Dilute	ELISA: User Optimized, WB: User Optimized
Application Notes	Surface Lipoprotein p27 is suitable as a control in immunological assays. Specific conditions for reactivity should be optimized by the end user. Expect a band at 73.3 kDa for p27-MBP, (30.9 kDa for p27 and 42.4 kDa for MBP) in size corresponding to Surf