



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Hemoglobin S Control Peptide ROC-000-001-GS5

Article Name	Hemoglobin S Control Peptide
Biozol Catalog Number	ROC-000-001-GS5
Supplier Catalog Number	000-001-GS5
Alternative Catalog Number	ROC-000-001-GS5
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Concentration	1.0 mg/ml by dry weight
Buffer	None
Form	Lyophilized
Application Notes	HbS peptide can be used as a control with the HbS antibody. Control peptide should be used at 1.0 $\mu$ g per 1.0 $\mu$ l of antiserum in per assay.