



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone H3, 21-35, Unconjugated ROC-000-001-K20

Article Name	Histone H3, 21-35, Unconjugated
Biozol Catalog Number	ROC-000-001-K20
Supplier Catalog Number	000-001-K20
Alternative Catalog Number	ROC-000-001-K20
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	Unconjugated
Concentration	1.0 mg/mL by dry weight
Buffer	None
Form	Lyophilized
Application Notes	Histone H3, 21-35 Control Peptide is suitable for use in ELISA, Western Blot, Dot blot, PCA, and other assays. Control peptide should be used at 1.0 μ g per 1.0 μ l of antiserum in per assay. Specific conditions for reactivity should be optimized by the en