



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

CEFT panel 23 peptides ROC-000-001-L23

Article Name	CEFT panel 23 peptides
Biozol Catalog Number	ROC-000-001-L23
Supplier Catalog Number	000-001-L23
Alternative Catalog Number	ROC-000-001-L23
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Concentration	1.0 mg/ml by dry weight
Buffer	None
Form	Lyophilized
Application Notes	T-Cell Assay Panel Control Peptides are suitable for use in ELISA, Western Blot, Dot blot, PCA, and other assays. Control peptide should be used at 1.0 μ g per 1.0 μ l of antiserum in per assay. Specific conditions for reactivity should be optimized by the