



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## CEFT panel 32 peptides ROC-000-001-L24

Article Name	CEFT panel 32 peptides
Biozol Catalog Number	ROC-000-001-L24
Supplier Catalog Number	000-001-L24
Alternative Catalog Number	ROC-000-001-L24
Manufacturer	Rockland Immunochemicals
Category	Kits/Assays
Concentration	1.0 mg/ml by UV absorbance at 280 nm
Buffer	None
Form	Lyophilized
Application Notes	T-Cell Assay Panel Control Peptides are suitable for use in ELISA, Western Blot, Dot blot, PCA, and other assays. Control peptide should be used at 1.0 $\mu$ g per 1.0 $\mu$ l of antiserum in per assay. Specific conditions for reactivity should be optimized by the