



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone H3, 21-44 Fluorescein conjugated, FITC ROC-000-002-K13

| Article Name | Histone H3, 21-44 Fluorescein conjugated, FITC |
|----------------------------|--|
| Biozol Catalog Number | ROC-000-002-K13 |
| Supplier Catalog Number | 000-002-K13 |
| Alternative Catalog Number | ROC-000-002-K13 |
| Manufacturer | Rockland Immunochemicals |
| Category | Proteine/Peptide |
| Conjugation | FITC |
| Concentration | 1.0 mg/mL by dry weight |
| Buffer | None |
| Form | Lyophilized |
| Application Notes | Histone H3 5CF, 21-44 Control Peptide is suitable for use in ELISA, Western Blot, Dot blot, PCA, and other assays. Control peptide should be used at 1.0 μ g per 1.0 μ l of antiserum in per assay. Specific conditions for reactivity should be optimized by th |