



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## HIV-1 tat, Fluorescein conjugated, FITC ROC-000-002-M47

Article Name	HIV-1 tat, Fluorescein conjugated, FITC
Biozol Catalog Number	ROC-000-002-M47
Supplier Catalog Number	000-002-M47
Alternative Catalog Number	ROC-000-002-M47
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Conjugation	FITC
Concentration	1.0 mg/mL by dry weight
Buffer	None
Form	Lyophilized
Application Notes	HIV-1 tat, 5FAM Control Peptide is suitable for use in ELISA, Western Blot, Dot blot, PCA, and other assays. Control peptide should be used at 1.0 $\mu$ g per 1.0 $\mu$ l of antiserum in per assay. Specific conditions for reactivity should be optimized by the end